
Dräger

EV 800

Elektronik-Ventilator Electronic Ventilator

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Contents

	Page		Page
For Your Safety and That of Your Patients	5	Maintenance Intervals	24
Intended Use	6	– Replacing Inlet Filter	24
Ventilation Monitoring	7	– Reading Operating Hours	24
Preparation	8–16	Care	25–29
– Assembly of Ventilation Hose System	8	– Stripping Down	25
– Fitting Heat and Moisture Exchanger	8	– Cleaning/Disinfection in Hospital	27
– Connecting Aquapor Humidifier	9	– Cleaning/Disinfection in the Home	28
– Connecting Spirolog 1 Volume Monitor	10	Test of Readiness for Operation	30–33
– Connecting PEEP Valve	11	– Checking Ventilation Function	30
– Power Supply	12	– Checking Mains Power Failure Alarm	30
– Checking Internal Battery Charge	12	– Checking Internal Battery	31
– Connecting to External 12 V DC Supply	13	– Checking End Expiratory Pressure	31
– Equipment Check	16	– Checking Trigger	31
Operation	17–23	– Checking Upper Pressure Alarm	32
– Controlled Ventilation IPPV-Assist	17	– Checking Lower Pressure Alarm	32
– Assisted Ventilation	18	– Leak Test	33
– Using PEEP	18	Fault - Cause - Remedy	34–35
– Pressure-controlled Ventilation	19	What's What	37–39
– Synchronised Intermittent Mandatory Ventilation SIMV	19	Technical Data	40–41
– Apnoea-Alarm	20	Functional Description of Pneumatic Function	42–43
– Using Increased O ₂ Concentration	20	Parts-List	44–47
– Mains Power Failure	22		
– Alarms	22		
– Shut-down	23		

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For Your Safety and That of Your Patients

For correct and effective use of the apparatus and to avoid hazards it is essential to read the following recommendations and to act accordingly¹⁾:

Strictly follow the instructions for use

Any use of the apparatus requires full understanding and strict observation of these instructions.

The apparatus is only to be used for purposes specified here.

Maintenance ²⁾

The apparatus must be inspected²⁾ and serviced²⁾ by experts at regular 6 month intervals (and a record kept).

We recommend obtaining a service contract with Dräger-Service.

Repairs²⁾ and general overhaul of the apparatus may only be carried out by DrägerService.

Only original Dräger spare parts may be used for maintenance.

Observe chapter "Maintenance Intervals".

Power connection

The requirements laid down in IEC 601/1 "Safety of Medical Electrical Equipment" are applicable for electrically powered equipment.

Not for use in areas of explosion hazard

This apparatus is neither approved nor certified for use in areas where combustible or explosive gas mixtures are likely to occur.

Connection with other electrical equipment

Electrical connections to equipment which is not listed in these instructions for use should only be made following consultations with the respective manufacturers or an expert.

Liability for proper function or damage

The liability for the proper function of the apparatus is irrevocably transferred to the owner or operator to the extent that the apparatus has been serviced or repaired by personnel not employed or authorized by DrägerService or when the apparatus was used in a manner not conforming to its intended use.

Drägerwerk Aktiengesellschaft cannot be held responsible for damage caused by non-compliance with the recommendations given above. The warranty and liability provisions of the terms of sale and delivery of Drägerwerk Aktiengesellschaft are likewise not modified by the recommendations given above.

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¹⁾

Insofar as reference is made to laws, regulations or standards, these are based on the legal system of the Federal Republic of Germany.

²⁾

In accordance with DIN 31 051:

<i>Inspection</i>	<i>= examination of actual condition</i>
<i>Service</i>	<i>= measures to maintain specified condition</i>
<i>Repair</i>	<i>= measures to restore specified condition</i>
<i>Maintenance</i>	<i>= inspection, service and, if applicable, repair</i>

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Time-cycled, volume constant ventilator for adults and for children weighing about 12 kg or more:


- controlled, assisted ventilation (IPPV-Assist)
- controlled, assisted ventilation (IPPV-Assist) with pressure control
- *synchronised, intermittent mandatory ventilation (SIMV).*

With monitoring for:

- airway pressure P_{aw}
- electrical supply

Do not use the machine in the presence of flammable gases or flammable anaesthetic agents.

The unit should only be used upon medical advice and in line with the parameter settings and monitoring procedures prescribed by a physician.

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EV 800 with Aquapor humidifier and trolley

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Appropriate ventilation monitoring (recommended by DGAI)*

Ventilation monitoring must be appropriate and must indicate undesirable changes in the following ventilation parameters:

- airway pressure
- inspiratory oxygen concentration
- expiratory breathing volume
- inspiratory breathing gas temperature, where applicable.

Changes in these parameters may be caused by:

- acute changes in the patient's condition
- incorrect settings and assembly
- faulty equipment
- failure of power or gas supply.

The EV 800 has a built-in monitor for airway pressure. For clinical use, and to conform to DIN 13254 "Ventilators", the following additional monitoring is required:

- expiratory minute volume (e.g. Spirolog 1)
- inspiratory oxygen concentration (e.g. Oxydig), when additional equipment for oxygen enrichment is used.

For ventilation during transport (or in hospital) where constant monitoring by qualified medical personnel can be guaranteed, such monitoring may not be necessary.

During home ventilation when the patient does not require continuous ventilation and is able to recognise insufficient ventilation himself, monitoring of the expiratory minute volume may not be necessary.

Safeguarding ventilation with independent manual ventilation facility

When the life-saving function of the ventilator can no longer be guaranteed because of a known fault, the patient must immediately be ventilated on independent equipment with PEEP and/or increased inspiratory O₂ concentration (e.g. Resutator 2000) where required.

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* DGAI = German Association for Anaesthesia and Intensive-Care Medicine

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


* Option

Preparation

- Place the EV 800 on a level surface, such as a bedside table, wall console or trolley.

Assembly of Ventilation Hose System

- 1 Insert sterile microbial filter 654 ST with the isoclic connection into the inspiratory port as far as it will go, until the disconnect mechanism engages.

For home ventilation, the one-way bacterial filter 84 09 716 may also be used.
- 2 Fit 1.5 m hose to bacterial filter and to expiratory valve.

- 3 Fit coloured hose for airway pressure measurement to the left port (black ring) and to port on the expiratory valve.
- 4 Fit transparent control hose to the right port (**white ring**) and to the **white port** inside the cap on the expiratory valve.

Fitting Heat and Moisture Exchanger*

For short-term use, such as transport ventilation of tracheotomy patients

- Fit heat and moisture exchanger to patient connection on the expiratory valve.
- Note time/date on heat and moisture exchanger.

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


Connecting Aquapor Humidifier*

as an alternative to a heat and moisture exchanger:

- for stationary long-term use

Prepare Aquapor using appropriate Instructions for Use.

- 1 Attach Aquapor bracket to the right hand rail and secure.
- 2 Fit double connector into Aquapor.
- 3 Fill Aquapor bowl to top mark with distilled water.

- 4 Insert sterile microbial filter 654 ST with the ISOclic connection into the inspiratory port as far as it will go, until the disconnect mechanism engages.

For home ventilation, the one-way bacterial filter 84 09 716 may also be used.

- Fit 0.6 m hose to bacterial filter and to inlet port on the Aquapor.
 - Connect two 0.6 m hoses to a water trap and fit one end to outlet port on the Aquapor.
 - 5 Fit adaptor for temperature sensor into the other end.
- Connect thermometer to adaptor.
- 6 Fit expiratory valve into the adaptor.

- 7 Fit coloured hose for airway pressure measurement to the left port (black ring) and to the port on the expiratory valve.
- 8 Fit transparent control hose to the right port (**white ring**) and to the **white port** in the cap on the expiratory valve.

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

Connecting Spirolog 1* Volume Monitor

A volume monitor must be used for stationary operation in a hospital.

- Place Spirolog 1 on a level, horizontal surface such as a wall console, or on a trolley with the EV 800.

Prepare Spirolog 1 using appropriate Instructions for Use.

- 1 Fit 30 mm/22 mm reducer to the outlet on the expiratory valve.
- Connect two 0.6 m hoses to a water trap.
- 2 Fit one end onto the reducing adaptor.
- 3 Fit flow sensor into the other end of the hose. Connect flow sensor and Spirolog 1 by cable.

- Fit the hose with the flow sensor into a clamp on the rail.

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

Connection of PEEP Valve*


Prepare PEEP valve

- Use 30 mm input connector.
- Insert axis of valve actuator into the middle hole of input connector and housing.
- Screw together input connector and housing.

- Fit PEEP valve with 30 mm cone to the outlet on the expiratory hose.

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When Spirolog 1 is used

- Use 22 mm input connector.
 - Screw together parts of PEEP valve.
- 1 Push rubber sleeve onto flow sensor.
 - 2 Push PEEP valve with 22 mm cone into rubber sleeve.

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
Power Supply

For stationary use connect to mains power supply.

For use in wheelchair or car, an external 12 V supply can be used.

The internal battery must not be used as a power source during transport. It is for use only in emergencies, i.e. mains power failure or during change-over from one supply source to another.

Connecting to mains power supply


- Voltage on selector switch must correspond to mains voltage.
- Mains plug must be connected to a properly installed earthed mains power socket. Socket installed in accordance with VDE 0100 or equivalent national regulations if in doubt installation should be checked by a trained serviceman.

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Checking Internal Battery Charge

The internal battery for emergency use is automatically recharged from the mains during normal operation.

Check the charge of the battery before using for the first time, and re-charge when required.

- Mains switch on the rear to »0«.
- 1 Operating switch to operating mode.
 - 2 Yellow light »Internal Battery« flashes.
 - 3 Keep »Battery Test« button pressed until indicator on the display ceases to move.
 - 4 Charge indicator: red area = not sufficient green area = sufficient

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

A 3 hour charge should give about 45 minutes operation.

The battery cannot be recharged from an external 12 V DC supply.

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Connecting to External 12 V DC Supply

For instance, for emergency power supply or for use in ambulances.

- Only use cable no. 84 09 722.
- Never place battery on top of the machine; place the battery as far away from inlet filter of the ventilator as possible.
- Connect the free ends of the cable to battery correctly: white or red to + black to -
- Fit connector at the back of the machine.

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Preparation

Checking External 12 V Supply

- Mains switch on the back to »0«.
- 1 Switch operating mode switch to an operating mode.
- 2 Yellow light »External Battery« is on.

If not lit:

- 3 Check whether 12 V circuit breaker has been tripped. Press button.
- 4 Keep »Battery Test« button pressed until indicator on the display ceases to move.
- 5 Charge indicator: green area = sufficient.

Or:

- Machine has switched to internal battery: external battery is flat.
- Charge the battery with a separate charger, following manufacturer's instructions.

The machine can be operated for deep cycle type about 24 hours with a fully charged 74 Ah battery which is in good condition, and for about 10 hours with a 34 Ah battery.


Older batteries, or batteries in poor condition will have shorter operating times even though the indicator is in the green area.

Check battery charge during operation:

- Press »Battery Test« button and observe indicator.

Have a spare supply available.

After use, recharge battery immediately with separate battery charger - do not recharge from EV 800!

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Equipment Check

– immediately before use.

Dräger EV 800 Serial No. _____	Check list for EV 800 • Familiarity with current Instructions for Use is absolutely essential • Delete if not applicable, complete where applicable	Date: _____ Signature _____
---	--	------------------------------------

Check before each use

WHAT	HOW	REQUIRED
Breathing system	Hose Humidifier Expiratory valve Flow sensor and PEEP valve Test lung	Fully and firmly connected
Electricial supply	Mains switch »I« Operating mode »IPPV-Assist« Mains switch »0« Mains switch »I«	All lights are on for about 2 s; double-tone alarm sounds, machine then ventilates. Green LED lit. Red alarm light »Power Switch Over« on, accompanied by single-tone alarm Machine ventilates: green LED lit.
Alarm »High Pressure«	Keep outlet on expiratory valve sealed Release outlet	Red alarm display and single-tone ceases Pressure drops immediately, alarm tone ceases
Alarm »Low Pressure/Apnea«	Remove test lung from expiratory valve	After 2 ventilation strokes: red alarm light »Low Pressure/Apnea« on
Leak-proof system	Seal patient connection with a plug Machine on Standby »⏻«. Keep »Alarm Silence/Reset« button and »Battery Test« Pressed and switch to »IPPV-Assist« Switch again to Standby »⏻«.	Red alarm lights come on one after another and go out again; green trigger light comes on when leak test successful. Green trigger display goes out.

Attention
only the upper 3 red
alarm lights come

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
Controlled Ventilation IPPV-Assist

For frequency = 6 ... 38 bpm.

- 1 Operating mode switch to
»IPPV-Assist«.


Pre-setting

- 2 Set »Volume« in accordance with
body weight in kg.
Recommended: 8 to 15 mL/kg.
- 3 »Breath Rate« = 12 bpm.
- 4 »Insp. Time« = 1.5 sec;
the resulting exp. time is never
lower than insp. time.
- 5 Turn »Trigger« to left stop.
- 6 »Low Pressure« initially to 2 mbar.
- 7 »High Pressure« initially to
60 mbar.

To prevent hypoxaemia add O₂ when
required.

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When patient is connected

- Adjust setting for the therapy and
the patient.
- 8 Determine maximum airway
pressure P_{aw}
 - 9 Set »Low Pressure« 5 mbar below
maximum airway pressure.
 - 10 Set »High Pressure« 10 mbar
above maximum airway pressure.

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Regular routines


- About every 2 hours remove con-
densed water from the ventilation
hoses.

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Assisted Ventilation



- 1 Set trigger pressure about 3 mbar **below** end expiratory pressure to prevent the equipment triggering itself.

When the patient commences inspiration:

- 2 Green light »Trigger« flashes; inspiration follows, controlled by the patient.

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Using PEEP

- Rotate cap of PEEP valve until the chosen scale value is level with the edge of the cap.
Illustrated example: approx. 3 mbar.
 - The PEEP indicated on the meter must correspond within ± 2 mbar of the PEEP set on the PEEP valve.
- 3 Set »Trigger Pressure« about 3 mbar lower than measured PEEP to prevent the machine triggering itself.
 - 4 Adjust »Low Pressure« and »High-Pressure«.

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
Pressure-controlled Ventilation

For frequency 6 ... 38 bpm.

Function as for »IPPV-Assist«. However, when airway pressure reaches the set upper limit value »High Pressure«, the machine ceases inspiration, independent of volume supplied.

Alarm »HighPressure« is triggered when set values are exceeded by approx. 10 mbar.


- 1 Operating mode switch to »Pressure Cycled«.
- 2 Set »High Pressure«.
- 3 Set »Volume« no more than 40% higher than required.
»Volume« =
lung compliance x P_{aw} max.
- 4 The pressure indicated must correspond with »High Pressure« setting. If there is a sudden pressure rise, the pressure peak may not be indicated.

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Synchronised Intermittent Mandatory Ventilation SIMV

For frequency 1 ... 38 bpm.

Combination of machine ventilation and spontaneous breathing, with spontaneous breathing in between the machine strokes. If the patient does not breathe spontaneously, the machine takes over ventilation.

- 5 Operating mode switch to »SIMV«.
- 6 Set »Frequency« lower than in other modes.
- 7 Set »Trigger Pressure«; adjust for use with PEEP.

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Apnoea Alarm

For frequency setting 1 to 5 bpm:

When airway pressure P_{aw} does not fall below the set trigger pressure within 20 seconds, ventilation will commence at a frequency of 6 bpm. The red alarm light will come on and the double-tone alarm will sound.

For frequency settings of 6 bpm and over:

If patient does not breathe spontaneously, there will be no alarm and the machine will continue to ventilate at the frequency set.

As soon as spontaneous breathing re-commences, the machine will continue to operate in the SIMV mode.

Using Increased O_2 Concentration

For O_2 concentration over 21% by vol.

When the patient is dependent on increased O_2 concentration, use an O_2 monitor.

Oxygen assists and intensifies combustion.

Do not smoke; no naked lights!


Keep O_2 fittings free from grease/oil.

- 1 Fit 30 mm/22 mm reducer to inlet filter on the back.
- 2 Hang O_2 reservoir from the rail.
- Fit 0.6 m hose to reducer and to reservoir.
- 3 Connect hose to a O_2 source with flow control valve and flow meter.

Setting O_2 concentration:

Set O_2 flow and minute volume to give the desired O_2 concentration (see chart on page 21).

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
Power Failure

In the event of mains power failure or DC power supply failure, the machine automatically switches to the internal battery.

The machine can then operate for about 30 to 60 minutes in favourable conditions.

- 1 The red alarm light »Power Switch Over« flashes, the single-tone alarm sounds.
- 2 Press »Alarm Silence/Reset«* button.
- 3 The yellow light »Internal Battery« flashes.
- A single tone sounds every 5 minutes.

When the battery will operate for only 15 more minutes:

- 4 The red alarm light »Low Battery« flashes, the double-tone alarm sounds.

Provide an adequate alternative power supply immediately.

Immediately after each use re-charge the internal battery:


- Connect machine to power supply.
- Mains switch to »I«, green light »AC/Battery Charge« is on.

The internal battery cannot be recharged from an external battery or an external charger.

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Alarms

- Note the alarm tones.
urgent = double tone
information = single tone
- Identify fault via flashing lights.
- Press »Alarm Silence/Reset« button: alarm tone is switched off for about 60 seconds.
- Rectify fault.

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Double-tone alarm and flashing alarm
displays:

urgent, take action immediately.

Possible causes:

- »Low Pressure/Apnea«
- »Low Battery«

Single-tone alarm and flashing alarm
lights:

take suitable action

possible causes:

- »High Pressure«
- »Setting Error«
- »Power Switch Over« = power failure
- mains or external 12 V supply failure

When a fault has been rectified:

- A single tone is automatically
switched off (except for power
failure).
- The alarm light continues to flash.
- A double tone continues.

To reset alarm light and double tone:

- Press »Alarm Silence/Reset«
button: the alarm is reset.

Shut-down

After disconnection of patient

- Operating mode switch to
Standby »  «.
- Mains switch on Aquapor to »0«.
- Switch off the additional oxygen
supply.

**For lengthy breaks in operation
and before transport**

- Mains switch on EV 800 to »0«.
- Disconnect all mains power plugs.
- Arrange mains power cable at the
rear of the appliance.
- Check battery charge every
3 months, and recharge when
required.

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Inlet filter	replace when soiled; when used out of doors check daily, otherwise check weekly.
Inspection and maintenance	every six months by trained servicemen

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
Clean and disinfect equipment and/or components prior to any maintenance procedure and before shipping it for repair!

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Replacing Inlet Filter


- Check weekly.
Replace if soiled.
When the machine is used for transport ventilation check daily.

A soiled filter may affect the performance of the machine.

- 1 Unscrew plastic cap and remove the outer grid.
- 2 Take out soiled filter.
- 3 The inner grid should remain in the machine.
- 2 Insert new filter.
- 1 Re-insert outer grid and screw plastic cap back on.

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

Reading Operating Hours

- 4 Press »Alarm Silence/Reset« button and »Battery Test« button simultaneously.
- 5 Multiply indicator reading by 100 = operating hours.

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
Stripping Down


- 1 Pull PEEP valve out of the rubber sleeve. Disconnect the plug on the flow sensor. Pull flow sensor out of hose and rubber sleeve.
- 2 Remove 20mm/22mm reducer from hose and from expiratory valve.
- 3 Detach control and measuring hoses from the ports on the machine and the expiratory valve.
- 4 Remove expiratory valve with temperature sensor adaptor from the hose.
- 5 Detach hose from Aquapor. Remove double connector from Aquapor.
- 6 Press the area marked "PRESS" on the isoclic connection of the microbial filter with thumb and index finger when removing filter from the port.
Do not remove isoclic connection without having pressed the area marked "PRESS" as this will lead to premature wear or damage to the locking mechanism.
- Detach hoses from water traps.
- 7 Remove O₂ reservoir from the rail and disconnect reducer from inlet filter.

Stripping down O₂ reservoir

- 8 Remove bag.
- 9 Remove hose.
- 10 Remove reducer from hose.

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
Stripping down Aquapor

- 1 Unscrew Aquapor bowl with float and empty the water.
- 2 Unscrew nut and remove lid.

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Stripping down expiratory valve

- Remove adaptor from temperature sensor.
- 3 Remove upper part.
 - 4 Screw out nut.
 - 5 Remove balloon from the upper part. Do not strip any further.

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Flow sensor

If contamination of the flow sensor is suspected:

- Dispose of flow sensor with household waste.

Otherwise:

- Re-use flow sensor as long as calibration can be carried out.

The sensor cannot be autoclaved and cannot be disinfected.

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Cleaning/Disinfection in Hospital

For disinfection use surface disinfectants.

To ensure material compatibility, disinfectants based on the following are recommended:

- aldehydes,
- alcohols,
- quaternary ammonium compounds.

To avoid damage to materials, disinfectants based on the following compounds should not be used:

- halogen-releasing compounds,
- strong organic acids,
- oxygen-releasing compounds.

For users in the Federal Republic of Germany we recommend only disinfectants contained in the list published by the DGHM*. The DGHM list contains details of the ingredients of disinfectants.

For countries not covered by the DGHM list, we recommend the preparations on the above list.

EV 800, Aquapor, Spirolog 1, Oxydig, cable, thermometer, O₂ connecting hose, O₂ reservoir

- Wipe disinfection, for instance with Buraton 10 F (Messrs. Schülke & Mayr, Norderstedt); follow manufacturer's instructions.

Hoses, expiratory valve, PEEP valve, patient part of Aquapor, water traps, reducer, double adaptor

- Bath disinfection with Helipur (Messrs. Braun, Melsungen); follow manufacturer's instructions.
- Rinse with clear water - preferably from a soft water source - and shake off water thoroughly.
- Dry parts **completely**. Water left in the control hose or in balloon of expiration valve may impair function.

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* DGHM = German Society for Hygiene and Microbiology

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Microbial filter 654 ST

Do not use bath disinfection methods!
Do not disinfect in Aseptor!
Do not sterilize with ethylene oxide!

- Wipe surface with a damp disposable cloth.
Sterilize filter in autoclave at a maximum of 134°C.
The filter may be sterilized 24 times.
Mark each sterilization in the appropriate area on the filter housing (e.g. with waterproof felt-tip pen).

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Cleaning/Disinfection in the home

EV 800, Aquapor, Spirolog 1, Oxydig, cable, O₂ connecting hose, O₂ reservoir

- Wipe off any visible dirt with a disposable cloth which has been soaked in detergent - liquids must not get into the machine.

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
Hoses, individual parts of expiration valve, individual parts of PEEP valve, patient part of Aquapor, thermometer, water traps, reducer, double adaptor

- Wash well in warm soapy water.
- Rinse well under running water, shake off water.
- Dry parts well.
then:
 - Immerse in a solution containing 1 part vinegar and 2 parts water and rinse well.
 - Rinse thoroughly under running water.
 - Shake off water to ensure that control and measuring hoses and balloon of expiratory valve are completely free of water. Water left in the control hose or in the balloon may affect the function.
- Dry all parts well.

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
After disinfection

Assemble expiratory valve:

- 1 Check balloon:
 - undamaged and evenly mounted
 - do not pull the actuator.
- Insert balloon in upper part and tighten with nut.
- Push upper part firmly onto T-piece
 - as far as it will go.

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Assemble PEEP valve:

- Insert axis of valve actuator into the middle hole of input connector and the housing.
- Screw together input connector and housing.

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Before re-use

- Assemble ventilation hose system, page 8.
- Check readiness of operation, page 30.

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

Tests of Readiness for Operation

To be carried out each time after cleaning.

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
- All lights are on for about 2 sec, the double-tone alarm sounds, the test lung is then filled and emptied alternating between inspiration and expiration.
- 8 The pressure indicator shows the airway pressure, alternating between inspiration and expiration.

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
Mains Power Failure Alarm

- Disconnect mains plug.
- 9 Red alarm light »Power Switch Over« flashes, single-tone alarm sounds.

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
Checking Internal Battery

Mains switch to »0«.

- 1 Yellow light »Internal Battery« flashes.
 - 2 Keep »Battery Test« button pressed until the indicator on the display ceases to move.
 - 3 Charge indicator: red area = not sufficient green area = sufficient, Otherwise re-charge, see page 13.
- Mains switch again to »I«.

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
Checking End Expiratory Pressure

- PEEP valve to 0 mbar:
- 4 Indicator shows (0 ± 2) mbar at the end of expiration.
- PEEP valve to 10 mbar:
- 4 Indicator shows (10 ± 2) mbar at the end of expiration.

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Checking Trigger

- PEEP valve to 0 mbar.
- 5 »Trigger Pressure« to about 3 mbar:
 - 6 green trigger lamp lights up regularly.
 - Inspiration is instantly triggered.
 - 5 Turn trigger pressure to left stop again.

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
Checking Alarm

»High Pressure«

- Keep outlet of expiratory valve sealed.

At 60 mbar:

- 1 Red alarm light »High Pressure« is on and single-tone alarm sounds.
- Open outlet of expiratory valve; the pressure drops immediately; the single-tone alarm ceases with the next ventilation stroke.
- 3 Press »Alarm Silence/Reset« button.
- 1 Alarm light »High Pressure« goes out.

Checking Alarm

»Low Pressure«

- Remove test lung from expiratory valve:

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Tests of Readiness for Operation

Leak Test

Carry out using the machine's self test.

- 1 Operating mode switch to Standby »  «.

- Seal patient part on Y-piece by using plug, for instance.

- 2 Keep both »Alarm Silence/Reset« and »Battery Test« buttons pressed **simultaneously** and:

- 1 switch operating mode switch to »IPPV-Assist«:

the self test commences.

- The red alarm lights come on one after another and go out when test is successful.

- 3 Green trigger light remains on when test is completed.

- 1 Operating mode switch again to Standby »  «; green trigger light goes out.

If at the end of the leak test red lights still flash:


- Test unsuccessful.

Light »Low Pressure/Apnea« flashing:
ventilation hose system not leak-tight.

Light »High Pressure« flashing:
safety valve has not opened.

If all tests were successful, the machine is ready for operation.

Test Spirolog 1 and O₂ monitor in accordance with appropriate Instructions for Use.

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Alarm	Cause	Remedy
»Low Pressure«	Hose system disconnected or not firmly connected Leak in hose system due to damage Low limit value setting below PEEP	Attach hose system firmly to the connection Replace damaged components Set 5 mbar below max. airway pressure
»Apnea«	With SIMV frequency lower than 6 bpm, airway pressure has not passed trigger pressure setting No spontaneous breathing Wrong setting	Select different operating mode Correct the setting
»Low Battery«	Internal battery low	Use external battery or connect to mains power supply; re-charge internal battery immediately (see page 13)
»High Pressure«	Blocked airways Kink in hose Upper limit value setting too low	Suction off Ensure that hoses are not kinked Adjust upper airway pressure setting 10 mbar higher
»Setting Error«	Incorrect setting, e.g. inspiration time greater than expiration time In operating mode »IPPV-Assist« and »Pressure Cycled« frequency range setting lower than 6 bpm »Low Pressure« setting above »High Pressure« Microprocessor fault	Correct the setting Operating mode switch to Standby »  «. Mains switch to »0« and then back to »I«. Select operating mode

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
Alarm	Cause	Remedy
»Power Switch Over«	Machine has switched from mains power to internal or external battery, or from external battery to internal.	Check that power supply is adequate. Recharge battery immediately (see page 13)
Continuous tone alarm, all red displays lit simultaneously	Machine has recognised micro-processor failure and switched to emergency ventilation, with frequency 15 bpm and inspiration time = 2 sec, expiration time = 2 sec, taking into account stroke volume setting	Use a spare machine when necessary and call DrägerService.

If the faults cannot be rectified:
use a spare machine and call DrägerService.

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

What's What

Front view

- 1 Red alarm lights
 - »Low Pressure/Apnea«
 - »Low Battery«
 - »High Pressure«
 - »Setting Error«
 - »Power Switch over«
- 2 »Alarm Silence/Reset« button
To suppress alarm sound for 60 sec.
To re-set alarm when fault is rectified.
To test alarm sound function:
- 3 Control panel
- 4 Cover to protect settings with sign to describe alarms
- 5 Carrying handle
- 6 Inspiration port, 22mm ISO-cone
- 7 Port for control hose of expiratory valve (white ring)
- 8 Port for airway pressure measuring hose (black ring)
- 9 Indicator for airway pressure and for operating hours, when
»Alarm Silence/Reset« and »Battery Test« buttons are pressed simultaneously.
Multiply display by 100 = operating hours
- 10 Scale to show battery charge, only operates when »Battery Test« button is pressed and mains switch at »0«
- 11 Green light to show trigger operation
- 12 Lights to show actual energy source
 - green lamp »AC/Battery Charge«
 - yellow lamp »External Battery«
 - flashing yellow lamp »Internal Battery« (combined with single tone every 5 min)
- 13 »Battery Test« button to show actual charge of internal or external battery

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
Back view

- 1 Rail for accessories
- 2 Reel for mains cable
- 3 Vent – do not seal
- 4 Holder for mains cable
- 5 Socket for external 12 V DC power source
- 6 Circuit breaker for external 12 V DC power source
- 7 Printer connector
- 8 Mains power switch with fuses
- 9 Selector switch for mains voltage
- 10 Safety valve – do not seal
- 11 Inlet filter – do not seal
- 12 Vent slits underneath – do not block

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Ambient conditions

During operation:	
Temperature	10 to 40°C
Air pressure	530 to 1100 hPa
Humidity	0 to 90% rel. humidity

During storage:	
Temperature	-20 to 70°C
Air pressure	530 to 1100 hPa
Humidity	0 to 90% rel. humidity

Performance Data

Operating parameters:

Control principle	Volume measuring time- cycled, volume-constant or pressure-controlled
-------------------	---

Ventilation frequency	6 to 20 bpm (SIMV 1 to 20 bpm) in steps of 1 bpm 22 to 38 bpm in steps of 2 bpm
-----------------------	---

Inspiration time	0.5 to 5.5 sec
------------------	----------------

Tidal volume	0.1 to 2.2 L
Tolerance	± 8% of set value, at least 30 mL

Minute volume in accordance with ISO	minimum: 1 L/min maximum: 28 L/min
---	---------------------------------------

Inspiration pressure limit	25 to 100 mbar
----------------------------	----------------

Trigger pressure	-10 to +10 mbar
------------------	-----------------

optional Positive end expiratory pressure PEEP	0 to 10 mbar
Tolerance	± 2 mbar

Machine compliance with 1.5 m hose	≤ 1 mL/mbar
with Aquapor humidifier and patient hose system	≤ 2.5 mL/mbar

Pressure indication	-10 mbar to 100 mbar
Tolerance	
Range -10 to +10 mbar	± 2 mbar
Range +10 to 100 mbar	± 5 mbar

Monitoring:

Airway pressure P_{aw}	
Low airway pressure	2 to 60 mbar
High airway pressure	25 to 100 mbar

Alarms	visual and audible audible alarm can be suppressed for 1 minute
--------	---

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Electrical power:

Drive	piston pump with brushless induction motor
Electrical connection	100 – 127 V 50/60 Hz with switch over to 220 – 240 V 50/60 Hz or 12 V DC voltage with external or internal battery
Connection for external 12 V direct voltage	minimum 11.6 V maximum 30 V
Current consumption	maximum 5 A
Operating time with internal battery	at least 45 minutes average 1 hour
Operating time with external battery, e.g. 74 Ah	at least 24 hours
Power input at 12 V DC voltage at 220 – 240 V / 100 – 127 V	50 W 110 W
Fuses	1 A fuse in mains switch 6 A circuit breaker for external 12 V supply
Safety class 1	 = DIN IEC 601, section 19, table 4: max. permissible earth leakage current 0.5 mA
Dimensions (W × H × D) mm	318 × 235 × 343
Weight	17 kg

Connection for printer

Centronics standard:


Pin	Description	Remark
1	DATA STROBE	
2	DATA 1	
.	.	
9	DATA 8	
11	BUSY	
12	PAPEROUT	
14	—	do not use
15	—	do not use
31	—	do not use
34	—	do not use

The connector is in direct electrical contact with the electronics of EV 800. The electronic chassis is at conductor potential.

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Key to pneumatic functional diagram

- 1 Filter
- 2 Non-return valve
- 3 Pump
- 4 Non-return valve
- 5 Safety valve
- 6 Non-return valve
- 7 3/2 way solenoid valve
- 8 Filter
- 9 Pressure transducer
- 10 Expiratory valve
- 11 PEEP valve
- 12 O₂ storage bag
- 13 Safety valve
- 14 Non-return valve

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Ventilation - Inspiration

Pump **3** supplies air via the non-return valve **4**, filter **8** and expiratory valve **10** to the patient.

Simultaneously, solenoid valve **7** controls and closes the expiratory valve **10**.

Airway pressure is measured by the pressure transducer **9**.

During the inspiration time, the expiratory valve **10** remains controlled, and the system is closed.

When airway pressure reaches a pre-set pressure limit, inspiration is immediately ended, and the pump stops. The solenoid valve **7** releases the control pressure on the expiratory valve **10**, which then opens - expiration can now take place.

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Functional Description of Pneumatic Function

Safety Valve

The safety valve **5** limits the maximum pressure in the machine to 100 mbar.

Ventilation - Expiration

The control pressure on the expiratory valve **10** is released by the solenoid valve **7**, the expiratory valve **10** opens and the patient can breathe out through this outlet.

Simultaneously, the pump **3** draws in ambient air via filter **1** and non-return valve **2**.

PEEP

PEEP valve **11** produces a positive end expiratory pressure.

Spontaneous breathing

During the expiratory phase, with the machine switched off or during power failure, spontaneous breathing is possible:

Inspiration takes place via filter **1**, non-return valve **6**, filter **8** and expiratory valve **10**.

Expiration takes place via the outlet on the expiratory valve **10**.

Additional air

In all operating modes, the patient can breathe in additional air via non-return valve **6** and filter **1**.

Added O₂

O₂ flows continuously into O₂ storage bag **12**.

The pump initially draws oxygen out of the storage bag; any additional volume required is drawn in via the non-return valve **14** from the ambient air. During O₂ failure, ventilation with ambient air takes place.

The safety valve **13** prevents the bag from overfilling.

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Parts List

* These parts are required in several numbers for regular servicing.

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No. in Fig. p. 47	Name and description	Code No.
	Basic unit EV 800 (100 – 240 V)	84 09 364
2	Bacterial filter (disposable)	84 09 716*
3	alternatively: Microbial filter 654 St (5 ea)	67 33 895*
4	Hose systems without Aquapor Hose system (disposable)	84 09 723*
5	alternatively: Patient/control hoses	84 09 719*
5.1	Hose, 1.5 m	84 09 894
5.2	Set of control hoses	84 09 891
6	Expiratory valve	84 09 690*
6.1	Mushroom valve	84 09 709
7	Optionally for item 5: PEEP Valve (Ambu PEEP)	84 07 475*
8	Set of catheter connectors E adult (12 ea)	84 03 685*
9	one each of size 6 to 12.5 Corrugated hose	84 02 041*
10	additionally required for volume measurement with item 5 Adapter 30/22 mm	84 09 717*
11	Hose, 1.5 m	84 09 894*
12	Hose system for Aquapor Hose system EV 800	
12.1	humidifier	84 09 982*
12.2	Hose, 0.6 m	84 09 718
12.3	Adapter 30/22 mm	84 09 717
12.4	Water trap	84 04 985
12.5	Set of control hoses	84 09 891
12.6	Sensor sleeve E (adult)	84 05 791
12.7	Thermometer	2M 13 259
13	Socket	M 25 647
13.1	Expiratory valve	84 09 690*
14	Mushroom valve	84 09 709
15	optionally for item 12 PEEP valve (Ambu PEEP)	84 07 475*
16	Set of catheter connectors E , adult (12 ea)	84 03 685*
17	one each of size 6 to 12.5 Corrugated hose	84 02 041*
17.1	Humidification Aquapor	84 05 020
17.2	Patient part for Aquapor	84 05 029*
18	Bowl	84 03 976
	Set of brackets	84 03 345

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
No. in Fig. p. 47	Name and description	Code No.
–	Equipment set	84 08 959
–	Aquapor holder	M 30 023
–	Monitor shelf	
–	Equipment set	M 30 058
–	depositing tray	M 30 060
–	Socket strip	
–	Miscellaneous accessories	
–	12 V cable with	84 09 722
–	angled connector	84 06 295
–	Hinged arm for rail	
–	Wearing parts	
–	Air filter EV 800	84 09 715

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These instructions for use apply
only to

EV 800

with serial-no.:

Without entry of serial-no. by Dräger
these instructions for use are pro-
vided for general information only and
are not destined to be used with a
specific device!

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GA 5665 d/e - 90 27 326

© Drägerwerk AG
2nd edition - May 1989
Subject to alteration